

# Gas Leak Detector

# LD239



Gas leaks can cause detector noise and baseline instability, shorten column life and waste expensive carrier gas. GL Sciences' LD239 leak detector is an easy way to quickly identify leaks in your system.

LD239 provides the most sensitive He leak detection available in an inexpensive, compact, easily charged (via USB), and easy to use configuration.

## Specifications

Detection Method	: Thermal conductivity
Target Gases	: Helium, CO <sub>2</sub> , Ar, Ne etc (gases that are largely different from referential gas) ※Not for flammable and corrosive gases
Sensitivity (He)	: [Standard range] 0.005mL/min at minimum [High range] 0.0005mL/min at minimum
Display	: LCD
Temperature Range	: 10 - 40 °C
Battery	: Lithium ion battery installed, rechargeable through a USB port, 2.5 hours for a full recharge. 5 hour battery life with continuous use
Dimension	: 50 (W) x 20 (D) x 111 (H) mm
Weights	: 95g
Accessories	: USB cable (1M), Sample gas filter (probe built-in), referential gas filter (built-in), instruction manual, warranty certificate

Description	Cat.No.
Gas Leak Detector LD239	2702-19340
Battery for LD239	2702-19341
USB cable	2702-19331
Sample Filter	2702-19333
Reference Filter	2702-19334

## Worldwide Ordering Information

### **GL Sciences, Inc. Japan**

22-1 Nishishinjuku 6-Chome  
Shinjuku-ku, Tokyo,  
163-1130, Japan  
Phone: +81-3-5323-6620  
Fax: +81-3-5323-6621  
Email: [world@glsc.co.jp](mailto:world@glsc.co.jp)  
Web: [www.glsciences.com](http://www.glsciences.com)

### **GL Sciences B.V.**

De Sleutel 9  
5652 AS Eindhoven  
The Netherlands  
Phone: +31 (0)40 254 95 31  
Fax: +31 (0)40 254 97 79  
Email: [info@glsciences.eu](mailto:info@glsciences.eu)  
Web: [www.glsciences.eu](http://www.glsciences.eu)

### **GL Sciences, Inc. USA**

4733 Torrance Blvd. Suite 255  
Torrance, CA 90503  
Phone: 310-265-4424  
Fax: 310-265-4425  
Email: [info@glsciencesinc.com](mailto:info@glsciencesinc.com)  
Web: [www.glsciencesinc.com](http://www.glsciencesinc.com)



### **International Distributors**

Visit our Website at  
[www.glsciences.com/distributors](http://www.glsciences.com/distributors)